

Quality time for better care



## Product Catalog (US)

Conventional Needles • Micro-Collection Products  
Safety Products • Winged Infusion Sets • I.V. Catheters

# Table of Contents

## Hypodermics

SurGuard™3 Safety Hypodermic Needles and Syringes	3-4
Terumo Conventional Hypodermic Needles	5
1cc Tuberculin Luer Slip Syringes	6
3cc, 5cc, 10cc Syringes with Needle	6
Syringes without Needle	7

## Micro-Collection

Surshield™ Safety Blood Collection Sets	8
Multi-Sample Luer Adapter	8
Capiject™ Safety Lancets	9

## Winged Infusion Needle Sets


Surshield™ Safety Winged Infusion Sets	10
SURFLO™ Winged Infusion Sets	10

## I.V. Catheters


SURFLO™ ETFE I.V. Catheters	11
SurFlash™ Polyurethane I.V. Catheters	12
Surshield-PUR Safety I.V. Catheters	13

## Frequently Asked Questions and Technical Support

14



# Pick how you click.




**With 3 ways to activate** and a broader range of sizes than most other brands, SurGuard™3 is the industry's most adaptable safety needle. Add separate thumb and finger grips, sharper needles, and competitive pricing, and you'll understand why we say SurGuard™3 is the one needle that fits *every* clinician's style.

- 3 modes of activation—thumb, finger, or hard surface
- Separate finger and thumb grips with ridges and nonslip surfaces
- Hub and needle locks for security post activation
- Broad range of sizes
- Firm injection feel
- Sharpest needles among major brands tested\*
  - 10% sharper than similar needles from Becton, Dickinson and Company
  - More comfortable injections for patients

**Bottom line:** No matter what you're looking for in a needle, you're looking for SurGuard™3.

\*Reference: Data on file.

# SurGuard™ 3 Safety Hypodermic Needles and Syringes



## SurGuard™ 3 Safety Hypodermic Needles

Packed: 100/box, 8 boxes/case (800)

Product Code	Description
SG3-1825	18G x 1" needle
SG3-1838	18G x 1 1/2" needle
SG3-2025	20G x 1" needle
SG3-2038	20G x 1 1/2" needle
SG3-2051	20G x 2" needle
SG3-2125	21G x 1" needle
SG3-2138	21G x 1 1/2" needle
SG3-2225	22G x 1" needle
SG3-2238	22G x 1 1/2" needle
SG3-2325	23G x 1" needle
SG3-2338	23G x 1 1/2" needle
SG3-2516	25G x 5/8" needle
SG3-2525	25G x 1" needle
SG3-2538	25G x 1 1/2" needle
SG3-2713	27G x 1/2" needle
SG3-3013	30G x 1/2" needle

## SurGuard™ 3 Hypodermic Syringes with Safety Needle

Packed: 100/box, 4 boxes/case (400)

Product Code	Description
SG3-01T2516	1CC Syringe with 25G x 5/8" needle
SG3-01T2713	1CC Syringe with 27G x 1/2" needle
SG3-03L2025	3CC Syringe with 20G x 1" needle
SG3-03L2038	3CC Syringe with 20G x 1 1/2" needle
SG3-03L2125	3CC Syringe with 21G x 1" needle
SG3-03L2138	3CC Syringe with 21G x 1 1/2" needle
SG3-03L2225	3CC Syringe with 22G x 1" needle
SG3-03L2238	3CC Syringe with 22G x 1 1/2" needle
SG3-03L2325	3CC Syringe with 23G x 1" needle
SG3-03L2516	3CC Syringe with 25G x 5/8" needle
SG3-03L2525	3CC Syringe with 25G x 1" needle

\*Data on file.

These products have no components made of natural rubber latex.



## Conventional Hypodermic Needles

Our double-bevel design is precision-ground and honed to exacting quality standards. This attention to detail enables Terumo needles to penetrate skin with less force, which can mean greater comfort for patients.

### Terumo Conventional Hypodermic Needles

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ NN1825R	18G x 1"
■ NN1838R	18G x 1 1/2"
■ NN1925R	19G x 1"
■ NN1938R	19G x 1 1/2"
■ NN2025R	20G x 1"
■ NN2038R	20G x 1 1/2"
■ NN2125R	21G x 1"
■ NN2138R	21G x 1 1/2"
■ NN2225R	22G x 1"
■ NN2238R	22G x 1 1/2"
■ NN2325R	23G x 1"
■ NN2338R	23G x 1 1/2"
■ NN2516R	25G x 5/8"
■ NN2525R	25G x 1"
■ NN2538R	25G x 1 1/2"
■ NN2613R	26G x 1/2"
■ NN2713R	27G x 1/2"
■ NN2732R	27G x 1 1/4"
■ NN3013R	30G x 1/2"
■ NN3025R	30G x 1"



\* In Independent laboratory testing, Terumo needles were shown, on average, to be the sharpest among brands tested:

10% sharper than BD (22G - 30G)  
16% sharper than Covidien (22G - 27G)

When you think hypodermics...  
Think sharp. Think Terumo.



\*Data on file.  
These products have no components made of natural rubber latex.

# Conventional Hypodermic Syringes with Needle



## 1cc Tuberculin Luer Slip Syringes

*Packed: 100/box, 10 boxes/case (1,000)*

Graduations: 0.01 cc (1 cc total)

Product Code	Description
SS-01T	1 cc TB syringe without needle, luer slip
SS-01T2516	1 cc TB syringe with 25G x 5/8" removable needle
SS-01T2609	1 cc TB syringe with 26G x 3/8" removable needle
SS-01T2713	1 cc TB syringe with 27G x 1/2" removable needle

## 3cc Syringes with Needle

*Packed: 100/box, 10 boxes/case (1,000)*

Style: Luer Lock Tip

Graduations: 0.1cc (3cc total)

Product Code	Description
SS-03L2025	3cc syringe with 20G x 1" needle
SS-03L2038	3cc syringe with 20G x 1 1/2" needle
SS-03L2125	3cc, syringe with 21G x 1" needle
SS-03L2138	3cc syringe with 21G x 1 1/2" needle
SS-03L2225	3cc syringe with 22G x 1" needle
SS-03L2238	3cc syringe with 22G x 1 1/2" needle
SS-03L2325	3cc syringe with 23G x 1" needle
SS-03L2338	3cc syringe with 23G x 1 1/2" needle
SS-03L2516	3cc syringe with 25G x 5/8" needle
SS-03L2525	3cc syringe with 25G x 1" needle



## 5cc Syringes with Needle

*Packed: 100/box, 6 boxes/case (600)*

Style: Luer Lock Tip

Graduations: 0.2cc (6cc total)

Product Code	Description
SS-05L2025	5cc syringe with 20G x 1" needle
SS-05L2038	5cc syringe with 20G x 1 1/2" needle
SS-05L2125	5cc syringe with 21G x 1" needle
SS-05L2138	5cc syringe with 21G x 1 1/2" needle
SS-05L2225	5cc syringe with 22G x 1" needle
SS-05L2238	5cc syringe with 22G x 1 1/2" needle

## 10cc Syringes with Needle

*Packed: 100/box, 6 boxes/case (600)*

Style: Luer Lock Tip

Graduations: 0.2cc (12cc total)

Product Code	Description
SS-10L2025	10cc syringe with 20G x 1" needle
SS-10L2038	10cc syringe with 20G x 1 1/2" needle
SS-10L2125	10cc syringe with 21G x 1" needle
SS-10L2138	10cc syringe with 21G x 1 1/2" needle

These products have no components made of natural rubber latex.

## Hypodermic Syringes

The plastic blister package is tamper evident for assurance of sterility. The peel-packs provide waste disposal cost savings versus the bulky plastic tube containers of some other syringes. Perhaps most importantly, Terumo syringes with needles offer maximum patient comfort delivered through ultimate sharpness and smoothness.



Luer Lock Tip



Eccentric Luer Tip



Catheter Tip

### 3cc Syringes

*Packed: 100/box, 10 boxes/case (1,000)*

*Graduations: 0.1cc (3cc total)*

Product Code	Description
Without Needles:	
SS-03L	Luer Lock Tip
SS-03S	Luer Slip Tip

### 20cc Syringes

*Packed: 50/box, 4 boxes/case (200)*

*Graduations: 1 cc (25cc total)*

Product Code	Description
Without Needles:	
SS-20L2	Luer Lock Tip
SS-20S2	Luer Slip Tip
SS-20ES2	20cc Eccentric Luer Tip

### 5cc Syringes

*Packed: 100/box, 6 boxes/case (600)*

*Graduations: 0.2cc (6cc total)*

Product Code	Description
Without Needles:	
SS-05L	Luer Lock Tip
SS-05S	Luer Slip Tip

### 30cc Syringes

*Packed: 25/box, 8 boxes/case (200)*

*Graduations: 1 cc (35cc total)*

Product Code	Description
Without Needles:	
SS-30L	Luer Lock Tip
SS-30S	Luer Slip Tip
SS-30ES	30cc Eccentric Luer Tip

### 10cc Syringes

*Packed: 100/box, 6 boxes/case (600)*

*Graduations: 0.2cc (12cc total)*

Product Code	Description
Without Needles:	
SS-10L	Luer Lock Tip
SS-10S	Luer Slip Tip
SS-10ES	10cc Eccentric Luer Tip

### 60cc Syringes


*Packed: 25/box, 4 boxes/case (100)*

*Graduations: 1cc (60cc total)*

Product Code	Description
Without Needles:	
SS-60L	Luer Lock Tip
SS-60ES	Eccentric Luer Tip
SS-60C	2 oz. Catheter Tip

These products have no components made of natural rubber latex.

# Surshield™ Safety Winged Blood Collection Sets



## 12" Tubing

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
MN*SVS21B30	21G x 3/4" needle
MN*SVS23B30	23G x 3/4" needle
MN*SVS25B30	25G x 3/4" needle



## 7" Tubing

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
MN*SVS21B18	21G x 3/4" needle
MN*SVS23B18	23G x 3/4" needle

# Multi-Sample Luer Adapter



Sterile and disposable, adapters fit into standard infusion sets transforming them into winged blood collection sets. This allows multiple draws into evacuated tubes. Accepts standard tube holders.

Packed: 100/box, 30 boxes/case (3,000)


Product Code	Description
XX*MN2000TE	VENOJECT™ Multi-Sample Luer Adapter (male).

These products have no components made of natural rubber latex.


Terumo's Surshield Safety Winged Needle Sets with one-handed activation provides a solution with a cost-effective, easy to use design. Press shield over the needle until it locks into place.

- Unique needle locking mechanism helps reduce the risk of needlestick injuries
- One-handed activation for safety and ease of use
- No change in needle insertion technique
- Simple, built-in safety design is easy to use

1. Flip safety shield toward needle. Remove the needle and apply digital pressure to the site.



2. With the wing and shield between your thumb and index finger, press the safety shield against a hard surface such as a bedside table until an audible click is heard.



3. Visually confirm activation of the safety feature.



Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.



## Capiject™ Safety Lancets

### High quality, single-use lancets offer the form, fit and function you've come to rely on

- **Never see or handle the sharp:**  
Automatic blade or needle retraction styles - the sharp edge is enclosed before and after use
- **No special staging required:**  
Ready to use right out of the box
- **One-step activation:**  
Convenient activation for precise and consistent results
- **Color-coded sizes:**  
Quick identification of the right product





*Packed: 200/box, 10 boxes/case (2,000)*

Product Code	Color	Description	Width/Gauge	Depth
200101	Pink	Blade	1.5mm	1.0mm
200102	Green	Blade	1.5mm	1.5mm
200103	Blue	Blade	1.5mm	2.0mm
200104	Orange	Needle	23G	2.25mm
200105	Purple	Needle	28G	1.25mm

These products have no components made of natural rubber latex.

# Winged Infusion Sets



## Surshield™ Safety Winged Infusion Sets

### 12" Tubing

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
SV*S19BL	19G x 3/4" needle
SV*S21BL	21G x 3/4" needle
SV*S23BL	23G x 3/4" needle
SV*S25BL	25G x 3/4" needle

### 3 1/2" Tubing

SV*S19BLS	19G x 3/4" needle
SV*S21BLS	21G x 3/4" needle
SV*S23BLS	23G x 3/4" needle
SV*S25BLS	25G x 3/4" needle

## SURFLO™ Winged Infusion Sets

### 12" Tubing

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
SV-18BLK	18G x 3/4" needle
SV-19BLK	19G x 3/4" needle
SV-21BLK	21G x 3/4" needle
SV-22BLK	22G x 3/4" needle
SV-23BLK	23G x 3/4" needle
SV-25BLK	25G x 3/4" needle

### 8" Tubing

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
SV*25EL	25G x 1/2" needle
SV*27EL	27G x 1/2" needle

### 3 1/2" Tubing

Packed: 100/box, 5 boxes/case (500)

Product Code	Description
SV*19BLS	19G x 3/4" needle
SV*21BLS	21G x 3/4" needle
SV*23BLS	23G x 3/4" needle
SV*25BLS	25G x 3/4" needle

These products have no components made of natural rubber latex.

## SURFLO™ ETFE I.V. Catheters

A sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue. Preview chamber in clear sure-grip hub gives immediate indication of proper placement. Flexible thinwall catheter with large inner diameter assures good blood flow. All this plus safe, easy-to-handle, tamper-evident packaging make SURFLO™ I.V. Catheters a great clinical value.



*Packed: 50/box, 4 boxes/case (200)*

Product Code	Description
■ SR-OX1451CA	14G x 2"
■ SR-OX1632CA	16G x 1 ¼"
■ SR-OX1651CA	16G x 2"
■ SR-OX1664CA	16G x 2 ½"
■ SR-OX1832CA	18G x 1 ¼"
■ SR-OX1851CA	18G x 2"
■ SR-OX1864CA	18G x 2 ½"
■ SR-OX2025CA	20G x 1"
■ SR-OX2032CA	20G x 1 ¼"
■ SR-OX2051CA	20G x 2"
■ SR-OX2225CA	22G x 1"
■ SR-OX2419CA	24G x ¾"


These products have no components made of natural rubber latex.

# SurFlash™ Polyurethane I.V. Catheters

With SurFlash™, visual confirmation of catheter placement occurs due to a unique flashback feature. Flashback appears between the specially grooved needle and the catheter, providing medical staff with a higher certainty of proper catheter placement. The specially grooved needle also improves the ability to visualize flashback in patients where venipuncture is more difficult, and because SurFlash™ is made of polyurethane, which is firm at room temperature and softer at body temperatures, patient comfort is enhanced without compromising proper technique.

*Packed: 50/box, 4 boxes/case (200)*

Product Code	Description
SR*FF1451	14G x 2"
SR*FF1651	16G x 2"
SR*FF1664	16G x 2 ½"
SR*FF1832	18G x 1 ¼"
SR*FF1851	18G x 2"
SR*FF1864	18G x 2 ½"
SR*FF2025	20G x 1"
SR*FF2032	20G x 1 ¼"
SR*FF2051	20G x 2"
SR*FF2225	22G x 1"
SR*FF2419	24G x ¾"



*Note:* The needle insertion angle in the illustration differs from actual insertion angle.

These products have no components made of natural rubber latex.

## Surshield-PUR Safety I.V. Catheters

Surshield-PUR Safety I.V. Catheters have a passive safety activation with SurFlash™, our exclusive catheter placement indicator. The passive safety device with automatic activation reduces risk of needlestick injury without adding extra steps to your procedures. SurFlash™ enhanced flashback feature provides clear visual confirmation of catheter placement, for increased confidence and first-stick success.\*

*Packed: 50/box, 4 boxes/case (200)*

Product Code	Description
■ SR*SFF1832A	18G x 1 ¼"
■ SR*SFF2025A	20G x 1"
■ SR*SFF2032A	20G x 1 ¼"
■ SR*SFF2225A	22G x 1"
■ SR*SFF2232A	22G x 1 ¼"
■ SR*SFF2419A	24G x ¾"



\*Data on file.  
These products have no components made of natural rubber latex.

# Frequently Asked Questions and Technical Support

## ■ How can I contact Terumo Medical Products?

Please visit our website, [www.TerumoTMP.com](http://www.TerumoTMP.com), for more information or phone Customer Care at 1-800-888-3786

## ■ How can I obtain samples?

Please visit our website, [www.TerumoTMP.com](http://www.TerumoTMP.com) to request samples.

## ■ How can I order Terumo products through a distributor?

To find a distributor, please see Terumo website: [www.TerumoTMP.com](http://www.TerumoTMP.com) and click on Customer Support to find a List of Authorized Terumo distributors. Distributors should call Customer Administration for all questions related to orders and invoices: (800) 888-3786

## ■ How can I obtain pricing?

You must check with your local distributor for pricing since our business is conducted through authorized distributors.

## ■ Where are you located?

The Terumo Medical Corporate office is located at:  
265 Davidson Avenue, Suite 320, Somerset, NJ 08873  
1-800-283-7866  
[www.TerumoTMP.com](http://www.TerumoTMP.com)

## ■ What is the shelf life of Terumo products?

Standard & Safety Syringes & Needles: 5 years  
Standard I.V. Catheters: 3 years  
Standard & Safety Blood Collection/Infusion: 3 years  
Capject Lancets: 5 years

## ■ What's the difference between SurFlash™ & SURFLO™ I.V. Catheters?

SurFlash is an I.V. Polyurethane Catheter which allows for a visual confirmation of catheter placement due to its unique flashback feature. SURFLO is our regular ETFE I.V. Catheter. Both are radio-opaque and these products have no components made of natural rubber latex.

[www.TerumoTMP.com](http://www.TerumoTMP.com)

**RX ONLY.** The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.



**TERUMO** and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURFLO, SURSHIELD, SURFLASH and VENOJECT are registered trademarks of Terumo Corporation, Tokyo, Japan and they are registered with the U.S. Patent & Trademark Office. CAPIJECT and SURGUARD are registered trademarks of Terumo Medical Corporation, Somerset, NJ. © 2024 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. Rev. 07/24 PM-08410

**TERUMO MEDICAL CORPORATION** 265 Davidson Avenue Somerset NJ 08873 U.S.A.