

Hypodermic Needle Fact Sheet



Terumo's SurGuard®3 safety hypodermic needle offers three modes of activation and a wide range of benefits. See how SG3 is combining safety, ease of use and simple standardization to make everyone comfortable.

EASE OF USE

- Versatile, intuitive safety mechanism that mimics 75-80% of safety hypo market*

COMPLIANCE

- One single intuitive device could ensure compliance with internal standards*

SAFETY

- New and improved hub and needle lock to help clinicians with naked sharps
- Safety mechanism minimizes ability to override

IMPROVED CLINICAL OUTCOMES

- T-Sharp™ Technology may reduce injection pain*


UNIVERSAL COLOR CODING GAUGE SIZE

18G	19G	20G	21G	22G
23G	25G	26G	27G	30G

Some of these sizes may not be available in all regions.

* Data on file.


TERUMO SHARPNESS SUPERIORITY VS COMPETITORS



† Hypodermic needles with and without vial penetration; allergy and insulin syringes with vial penetration.

STANDARDIZATION IS KEY


Products from different hypodermic manufacturers can be consolidated into one standard, easy-to-use product.



VERSATILITY

Most syringes can be activated two ways but **SurGuard®3** has

3



Thumb, finger or hard surface

SurGuard®3's double-bevel design allows the needle to penetrate skin with less force, which could mean **greater comfort to patients***.



SurGuard®3
SAFETY HYPODERMIC NEEDLE

62.2%
of needle stick injuries
are due to discarded unprotected sharps*

RX ONLY. The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

TERUMO, TERUMO and SurGuard are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. PM-06400